## THE SMART AND PORTABLE LAB THAT IS -**CHANGING WINE ANALYSIS** IN YOUR CELLAR



The analysis platform developed to meet the needs of wineries and winemakers









## DNAPhone S.r.l.

Legal address: Viale Mentana 150 - 43121 Parma

Operational Headquarters: Via Ravasini 7 - 43126 Parma Tel. +39 0521 062397 - info@dnaphone.it

www.dnaphone.it









## APP AND CLOUD

PATENTED FIBER OPTIC
TECHNOLOGY AND
LED SOURCES

The extremely **intuitive** APP is used from the tablet and **guides the user throughout the entire analysis.** 

The system is automatically updated through the APP.


The **cloud** service allows the control and management of analysis results **anywhere**, **anytime**.

The **remote diagnostic** service ensures constant monitoring of the instrumentation in real time.



**Smart Analysis** is a portable, solid and compact device.

- •Dimension: **85X75X200 mm**
- •Weight: **1 kg**





wi-fi and bluetooth



thermostated up to **40° C** 



internal battery **7 Ah** + **microUSB** 



no **maintenance** 



spectrum **340-760 nm** 



The **simplified** analysis **methods** minimize the required manual operations.

The reagents are **pre-infused in disposable cuvettes**, associated with a dropper dispenser.



PARAMETER		RANGE
D-Glucose + D-Fruttose	0.10-8.00* g/l *extended to 320 g/l with "dilution kit"	
L-Lactic Acid	0.10-3.00 g/l	
L-Malic Acid		0.10-4.50 g/l
Acetic Acid		0.05-1.00 g/l
Total Acidity		0.50-8.00 g/l
Total SO <sub>2</sub>		20-300 mg/l
Free SO <sub>2</sub>		5-50 mg/l
Ammoniacal Nitrogen		20-140 mg/l
Alpha-amino Nitrogen		20-300 mg/l
Tartaric Acid		0.50-6.00 g/l
На		2.80-4.01
Anthocyanins		-
Ethanol	UP TO 15 % VOL /VOL *in combination with "dilution kit"	
Color	• Reading at 420-520-620 nm (OIV-MA-AS2-07B) • CIELab color space (OIV-MA-AS2-11)	

WINE AND GRAPE MUST ANALYSIS

The QR code allows **automatic system** calibration and data traceability.

